



Nieka[®] G-Series

The solution for high-throughput fusion laboratories

PERFECT FUSION EQUALS PREMIUM RESULTS

A superior control of the fusion process is the starting point
of first-class analytical results



AT NIEKA®, we used our strong experience from the fusion industry and forged an instrument that is the perfect fit for the most demanding high-throughput fusion laboratories. Downtime is your enemy and using the best tool is your only way to remain on top and deliver results on time.

For the last years, we spent thousands of hours creating just what you need: a high-throughput instrument that creates perfectly homogenous samples for analysis with no surprise. Thanks to the Nieka® exclusive monitoring features, you avoid the unpredictable, even if you are off-site.

For the first time in a gas fluxer, the G-Series offers a reliable temperature monitoring of the fusion process and fully independent fusion chambers (*patented*). Nothing has been left aside, from perfect monitoring to perfect fusion: this will lead you to premium analytical results.

HIGH PRODUCTIVITY:

Produces 4 to 8 glass disks or solutions per cycle (up to 48 samples per hour).

AUTOMATED OPERATION:

Multilingual touchscreen controller with easy to access controls, import and export fusion method from USB or network, easy to follow fusion program execution through the instrument's display or PC. All fusion parameters can be fully customized for optimal sample preparation.

MONITORING:

Temperature monitoring of the fusion process (*patent pending*), real-time monitoring through the Nieka® data management system.

EFFICIENCY AND REPEATABILITY:

Proprietary heating chambers (*patent pending*) allow a high thermal efficiency and higher fusion temperatures for difficult samples. Each chamber is well separated from the surroundings, allowing a perfect repeatability.

EASY INSTALLATION:

Only a low-pressure gas line and a standard electrical outlet are required. The G-Series runs on propane, any other liquid petroleum gas (LPG) or natural gas. No compressed air or oxygen is required.

OPTIONS:

The G-Series fluxers are available in 2 versions: 4 and 8 positions. Each unit can be delivered to make glass disks or solutions only or both (hybrid).

SAFETY:

Integrated locking safety door, real-time flame monitoring and automatic shut-off features, regulatory shutdown button.

CONNECTIVITY:

USB ports, TCP/IP Ethernet port.

REQUIREMENTS:

100-240V (50/60Hz) electrical supply, propane/ LPG (37 mbar - 0.5psi) or natural gas (20 mbar - 0.3 psi) line.



Built-in safety cabinet

You can see the entire fusion process and a locked door, that opens only when the cooling step is over, protects you. From the real-time flame monitoring to ultimate user safety, nothing has been left aside to provide you with the safest operation.

Sturdy mechanical systems

Nothing lasts forever but we're getting pretty close with the G-Series: every component is oversized and built to outlast anything else in the industry.

Fusion temperature measurement

Our proprietary fusion enclosure allows the instrument to run on low energy while achieving the required high fusion temperature for difficult samples. For the first time on a gas fluxer, the G-Series reports an accurate temperature measurement of the fusion process (*patent pending*).

Huge built-in touchscreen

Use the built-in touchscreen display to edit your methods easily and follow the fusion process in real-time.



Connectivity

USB or Internet, your fusion recipe can be uploaded to your instrument in a breeze. Manage your methods using our customized computer and tablet interface.

Advanced self-monitoring

Thanks to its self-monitoring feature, the G-Series fluxer analyzes all instrument functionality and logs all the cycles so we can handle the maintenance schedule remotely. That's one more thing off your shoulders.

Dedicated mold heating

G-Series fluxers have separate burners for the molds. You can pre-heat any size of mold while reducing your overall fusion process time. The hybrid includes a snap-in module to prepare solutions for AA or ICP sample preparation.

The G1: a basic line for every laboratory

NEW PRODUCT
available oct. 2018



SMALL AND EFFECTIVE

Produces up to 6 glass disk
per hour.

PACKED WITH ALL THE ESSENTIALS:

Same stirring technology and touchscreen controller as the GS instruments with 4 temperature levels and fully configurable agitation and cooling. The instrument also features our exclusive heating fusion chamber.

Nieka

4715 des Replats, Suite 135
Quebec, Canada, G2J 1B8
info@nieka.systems
www.nieka.systems

Tel. +1 418 476 1923